

Product Datasheet

Human Apelin 12 (AP12) ELISA Kit (orb1667868)

Description	Human Apelin 12 (AP12) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 93.7-6000 pg/mL. Sensitivity: 25.16 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.