

## Product Datasheet

# Human Angiotensin III (AngIII) ELISA Kit (orb1667931)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Angiotensin III (AngIII) ELISA Kit  |
| <b>Reactivity</b>          | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: 7.81-500 pg/mL.<br>Sensitivity: 2.77 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.<br>Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.   |