

Product Datasheet

Human Angiotensin II (AngII) ELISA Kit (orb1667932)

Description Human Angiotensin II (AngII) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 3.12-200 pg/mL.
Sensitivity: 1.13 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.

Expiration Date Please enquire.