

Product Datasheet

IKKB (phospho-S733) antibody (orb166795)

Description	Rabbit polyclonal antibody to IKKB (phospho-S733).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	IKBKB
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	86564
Uniprot ID	O14920
NCBI	NP_001177649.1 , NP_001229707.1 , NP_001547.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

