

Product Datasheet

Canine Albumin (ALB) ELISA Kit (orb1667995)

Description	Canine Albumin (ALB) ELISA Kit
Reactivity	Canine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.56-100 µg/mL. Sensitivity: 0.49 µg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.