

Product Datasheet

Bovine Alpha-2-Heremans Schmid Glycoprotein (aHSG) ELISA Kit (orb1668010)

Description	Bovine Alpha-2-Heremans Schmid Glycoprotein (aHSG) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.5-800 ng/mL. Sensitivity: 4.9 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)