

Product Datasheet

General Androsterone (ADT) ELISA Kit (orb1668042)

Description	General Androsterone (ADT) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 1.4-90 ng/mL. Sensitivity: 0.507 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.