

Product Datasheet

Rat Adrenocorticotrophic Hormone Receptor (ACTHR) ELISA Kit (orb1668116)

Description	Rat Adrenocorticotrophic Hormone Receptor (ACTHR) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.234 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates or other biological fluids.
Expiration Date	Please enquire.