

Product Datasheet

Human Angiotensin I Converting Enzyme 2 (ACE2) ELISA Kit (orb1668146)

Description	Human Angiotensin I Converting Enzyme 2 (ACE2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 7.7 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.