

Product Datasheet

Human ATP Binding Cassette Transporter C2 (ABCC2) ELISA Kit (orb1668164)

Description	Human ATP Binding Cassette Transporter C2 (ABCC2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.067 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 μ L. Sample Type: tissue homogenates, cell lysates or other biological fluids.
Expiration Date	Please enquire.