

Product Datasheet

Human Amyloid Beta Peptide 1-40 (Ab1-40) ELISA Kit (orb1668178)

Description	Human Amyloid Beta Peptide 1-40 (Ab1-40) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 5.58 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.