

Product Datasheet

Human Alpha-2-Plasmin Inhibitor (a2PI) ELISA Kit (orb1668185)

Description	Human Alpha-2-Plasmin Inhibitor (a2PI) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 12.5-800 ng/mL. Sensitivity: 4.3 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 100 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.