

Product Datasheet

TSC1 (phospho-Y304) antibody (orb166819)

Description	Rabbit polyclonal antibody to TSC1 (phospho-Y304).
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	Tsc1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	128746
Uniprot ID	Q9EP53
NCBI	NP_075025.2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

