

Product Datasheet

General 5-Hydroxytryptamine (5-HT) ELISA Kit (orb1668193)

Description	General 5-Hydroxytryptamine (5-HT) ELISA Kit
Reactivity	All
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 15.6-1000 ng/mL. Sensitivity: 4.16 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma or other biological fluids.
Expiration Date	Please enquire.