

Product Datasheet

TSC2 (phospho-S1096) antibody (orb166820)

| | |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Rabbit polyclonal antibody to TSC2 (phospho-S1096). |
| Species/Host | Rabbit |
| Reactivity | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | DOT |
| Immunogen | Synthetic Peptide |
| Target | Tsc2 |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 202071 |
| Uniprot ID | Q61037 |
| Dilution Range | DOT: 1:500 |
| Expiration Date | 12 months from date of receipt. |

