

Product Datasheet

Rat Phosphorylated microtubule-associated protein tau (pMAPT /pTAU) ELISA Kit (orb1668305)

Description	Rat Phosphorylated microtubule-associated protein tau (pMAPT /pTAU) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.066 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.