

Product Datasheet

ATG13 (phospho-S355) antibody (orb166833)

Description	Rabbit polyclonal antibody to ATG13 (phospho-S355).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	ATG13
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56572
Uniprot ID	O75143
NCBI	NP_001192048.1 , NP_001192051.1 , NP_001192050.1 , NP_001136145.1 , NP_055556.2 , NP_001192049.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

