

Product Datasheet

LC3 (phospho-S12) antibody (orb166845)

Description	Rabbit polyclonal antibody to LC3 (phospho-S12)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	This MAP1LC3A Antibody is generated from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding S12 of human MAP1LC3A.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q9GZQ8 , Q9H492
NCBI	NP_852610.1 , NP_115903.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

