

Product Datasheet

Human Receptor II for the Fc region of immunoglobulin G (FcGR2) ELISA Kit (orb1668454)

Description	Human Receptor II for the Fc region of immunoglobulin G (FcGR2) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 46.88-3000 pg/mL. Sensitivity: 17.35 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.