

Product Datasheet

NFkB p65 (phospho-S536) antibody (orb166852)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to NFkB p65 (phospho-S536) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | DOT |
| Immunogen | Synthetic Peptide |
| Target | RELA |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 60219 |
| Uniprot ID | Q04206 |
| NCBI | NP_001138610.1 , NP_068810.3 , NP_001230913.1 |
| Dilution Range | DOT: 1:500 |
| Expiration Date | 12 months from date of receipt. |

