

Product Datasheet

Bovine Lipopolysaccharide binding protein (LBP) ELISA Kit (orb1668560)

Description	Bovine Lipopolysaccharide binding protein (LBP) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.288 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.