

Product Datasheet

Bovine Neuropeptide Y (NPY) ELISA Kit (orb1668562)

Description	Bovine Neuropeptide Y (NPY) ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 24.69-2000 pg/mL. Sensitivity: 9.35 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.