

Product Datasheet

Rat Apolipoprotein C-II (APOC2) ELISA Kit (orb1668590)

Description Rat Apolipoprotein C-II (APOC2) ELISA Kit

Reactivity Rat

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.312-20 ng/mL.
Sensitivity: 0.137 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L.
Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.