

Product Datasheet

CCNB1 (phospho-S35) antibody (orb166866)

Description Rabbit polyclonal antibody to CCNB1 (phospho-S35).

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications DOT, WB

Immunogen Synthetic Peptide

Target CCNB1

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 48337

Uniprot ID [P14635](#)

NCBI [NP_114172.1](#)

Dilution Range WB: 1:1100, DOT: 1:500

Expiration Date 12 months from date of receipt.



WB analysis
of HeLa
cells(from
left to r...