

Product Datasheet

CCND3 (phospho-S264) antibody (orb166868)

Description	Rabbit polyclonal antibody to CCND3 (phospho-S264).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	CCND3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	32520
Uniprot ID	P30281
NCBI	NP_001129597.1 , NP_001129489.1 , NP_001751.1 , NP_001129598.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

