

Product Datasheet

Bird Influenza A virus subtype H5N1 IgG antibody (AIV-H5N1Ab-IgG) ELISA Kit (orb1668697)

Description	Bird Influenza A virus subtype H5N1 IgG antibody (AIV-H5N1Ab-IgG) ELISA Kit
Reactivity	Aves
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.