

Product Datasheet

Rat Anti-gonadotropin-releasing hormone receptor antibody (GNRHR-Ab) ELISA Kit (orb1668739)

Description	Rat Anti-gonadotropin-releasing hormone receptor antibody (GNRHR-Ab) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: Request Information. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.