

Product Datasheet

TSC1 (phospho-Y312) antibody (orb166877)

Description	Rabbit polyclonal antibody to TSC1 (phospho-Y312).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	TSC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	129767
Uniprot ID	Q92574
NCBI	NP_001155898.1 , NP_000359.1 , NP_001155899.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

