

## Product Datasheet

# Human Anti-rubella virus IgG antibody (RV-Ab-IgG) ELISA Kit (orb1668862)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Anti-rubella virus IgG antibody (RV-Ab-IgG) ELISA Kit   |
| <b>Reactivity</b>          | Human   |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: Detection antibody. Sensitivity: Request Information. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.   |