

Product Datasheet

Mouse Alpha-Melanocyte Stimulating Hormone (A-MSH) ELISA Kit (orb1668900)

Description	Mouse Alpha-Melanocyte Stimulating Hormone (A-MSH) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 123.5-10000 pg/mL. Sensitivity: 43.9 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.