

Product Datasheet

Rat Apelin 13 (AP13) ELISA Kit (orb1668919)

Description	Rat Apelin 13 (AP13) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 74.07-6000 pg/mL. Sensitivity: 26.38 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.