

Product Datasheet

Rat Fibroblast Growth Factor 10 (FGF10) ELISA Kit (orb1669003)

Description	Rat Fibroblast Growth Factor 10 (FGF10) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.81-500 pg/mL. Sensitivity: 3.3 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.