

Product Datasheet

Human Free Thyroxine (fT4) ELISA Kit (orb1669011)

Description	Human Free Thyroxine (fT4) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 4-64 pmol/L. Sensitivity: 2.2 pmol/L. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.