



Product Datasheet

ENO1 antibody (orb166904)

Description	Mouse monoclonal antibody to ENO1.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	ENO1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG (k)
Clonality	Monoclonal
MW	47169
Uniprot ID	P06733
NCBI	NP_001419.1 , NP_001188412.1
Dilution Range	IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...