




Product Datasheet

ERK1/2 antibody (orb166906)

Description	Mouse monoclonal antibody to ERK1/2
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Recombinant Protein
Target	MAPK3
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	43136
Uniprot ID	P27361
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis of extracts from HCC cell li...



WB analysis of (lane 1) Jurkat whole cel...