




Product Datasheet

VEGFR3 antibody (orb166909)

Description	Mouse monoclonal antibody to VEGFR3.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant Protein
Target	FLT4
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	152757
Uniprot ID	P35916
NCBI	NP_891555.2 , NP_002011.2
Dilution Range	IHC-P: 1:25, FACS: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Immunohistochemical analysis of formalin...