

## Product Datasheet

# Human Undercarboxylated Osteocalcin (ucOC) ELISA Kit (orb1669278)

<b>Description</b>	Human Undercarboxylated Osteocalcin (ucOC) ELISA Kit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.625-40 ng/mL. Sensitivity: 0.252 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids.
<b>Expiration Date</b>	Please enquire.