

Product Datasheet

PLM (phospho-S68) antibody (orb166931)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to PLM (phospho-S68). |
| Species/Host | Mouse |
| Reactivity | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Synthetic Peptide |
| Target | Fxyd1 |
| Form/Appearance | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS. |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles |
| Note | For research use only |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| MW | 10323 |
| Uniprot ID | Q9Z239 |
| NCBI | NP_062376.2 , NP_443717.1 , NP_919302.1 , NP_443718.1 |
| Dilution Range | WB: 1:100-500 |
| Expiration Date | 12 months from date of receipt. |

