

Product Datasheet

Human Soluble cluster of differentiation 14 (sCD14) ELISA Kit (orb1669324)

Description	Human Soluble cluster of differentiation 14 (sCD14) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 3.75-60 ng/mL. Sensitivity: 1.52 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 µL. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.