

Product Datasheet

Rat Glycine receptor subunit alpha-2 (GLRA2) ELISA Kit (orb1669384)

| | |
|----------------------------|---|
| Description | Rat Glycine receptor subunit alpha-2 (GLRA2) ELISA Kit |
| Reactivity | Rat |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 0.312-20 ng/mL. Sensitivity: 0.115 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids. |
| Expiration Date | Please enquire. |