

Product Datasheet

Mouse Interleukin 17 (IL-17) ELISA Kit (orb1669394)

Description	Mouse Interleukin 17 (IL-17) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 7.8-500 pg/mL. Sensitivity: 3.0 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.