

Product Datasheet

Human Interleukin 1 β (IL-1 β) ELISA Kit (orb1669397)

Description Human Interleukin 1 β (IL-1 β) ELISA Kit

Reactivity Human

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 15.6-1000 pg/mL.
Sensitivity: 5.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.