



Product Datasheet


PRDM2 (RIZ1) antibody (orb166940)

Description	Mouse monoclonal antibody to PRDM2 (RIZ1).
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Target	PRDM2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	188915
Uniprot ID	Q13029
NCBI	NP_056950.2 , NP_001129082.1 , NP_036363.2 , NP_001007258.1
Dilution Range	WB: 1:1000, ChIP: 1:100
Expiration Date	12 months from date of receipt.



ChIP analysis of MCF7 cells using PRDM2

...



ChIP analysis of MCF7 cells using PRDM2

...