

Product Datasheet

Human NACHT, LRR and PYD domains-containing protein 1 (NLRP1) ELISA Kit (orb1669435)

Description	Human NACHT, LRR and PYD domains-containing protein 1 (NLRP1) ELISA Kit
Reactivity	Human
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 78.1-5000 pg/mL. Sensitivity: 40 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.
Expiration Date	Please enquire.