



## Product Datasheet

### UBC9 antibody (orb166944)

<b>Description</b>	Mouse monoclonal antibody to UBC9.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Target</b>	UBE2I
<b>Form/Appearance</b>	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	18007
<b>Uniprot ID</b>	<a href="#">P63279</a>
<b>NCBI</b>	<a href="#">NP_919237.1</a> , <a href="#">NP_003336.1</a> , <a href="#">NP_919236.1</a> , <a href="#">NP_919235.1</a>
<b>Dilution Range</b>	WB: 1:100-500
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of mouse spleen ti...