

## Product Datasheet

# Rat Prostaglandin E2 (PG-E2) ELISA Kit (orb1669453)

**Description** Rat Prostaglandin E2 (PG-E2) ELISA Kit

**Reactivity** Rat

**Tested Applications** ELISA

**Note** For research use only

**Application notes** Detect Range: 3.12-200 pg/mL.  
Sensitivity: 1.04 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.  
Sample Volume: 1-200  $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids.

**Expiration Date** Please enquire.