

Product Datasheet

Bovine Progesterone (PROG) ELISA Kit (orb1669470)

Description Bovine Progesterone (PROG) ELISA Kit

Reactivity Bovine

Tested Applications ELISA

Note For research use only

Application notes Detect Range: 0.78-50 ng/mL.
Sensitivity: 0.4 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm.
Sample Volume: 1-200 μ L. Sample Type: Serum, Plasma, Other biological fluids.

Expiration Date Please enquire.