




Product Datasheet

MMP14 antibody (orb166952)

Description	Mouse monoclonal antibody to MMP14 (C-term).
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	MMP14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgM
Clonality	Monoclonal
MW	65894
Uniprot ID	P50281
NCBI	NP_004986.1
Dilution Range	WB: 1:100-500, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



WB analysis of A2058 whole cell lysates ...