

## Product Datasheet

# Chicken Corticosterone (CORT) ELISA Kit (orb1669579)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Chicken Corticosterone (CORT) ELISA Kit   |
| <b>Reactivity</b>          | Gallus  |
| <b>Tested Applications</b> | ELISA   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Detect Range: 0.05-10 ng/mL.<br>Sensitivity: 0.025 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 1-200 $\mu$ L. Sample Type: Serum, Plasma, Other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.   |