

Product Datasheet

TNK1 antibody (orb166967)

Description	Mouse monoclonal antibody to TNK1
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	TNK1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	72468
Uniprot ID	Q13470
NCBI	NP_001238831.1 , NP_003976.2
Dilution Range	WB: 1:100-1000
Expiration Date	12 months from date of receipt.

